



Services and Techniques

Eurofins EAG Laboratories in Toulouse, France, is an independent third-party laboratory. We provide our clients with local, scientifically sound solutions and fast turnaround times, including 72 to 24 hours upon receipt of the sample. Your satisfaction is our priority.

New: Glow Discharge Optical Emission Spectrometry (GDOES)

Surface Analyses and Depth Profiling

Using a customer-focused development approach, we added a new technique to our characterization services – Glow Discharge Optical Emission Spectroscopy, or GDOES. This technology relies on the application of Argon plasma on the surface of a material to perform sputtering. The capacity of the plasma to be pulsed, while being controlled by a Radio Frequency generator, unlocks the possibility of characterization of thin layers, with enhanced depth resolution thanks to slower sputtering rates.

Analyze Wide Spectrum of Wavelengths Simultaneously

The analysis itself relies on the acquisition of the intensity of light obtained from deexcitation of the atoms and ions in the plasma. A polychromator is used to




Figure 1: Basic scheme extracted from Horiba GD-Profilier 2 brochure

analyze a wide spectrum of wavelengths simultaneously and strengthens the use of the pulsed plasma with the goal of reaching very good depth resolutions, down to nanometer depending on uses and analytical needs.




Figure 2 : Observation of a crater obtained after GDOES sputtering

25 Years of Expertise

This technological addition uses our strong knowledge of plasma behavior and completes our instrumental capacity on elemental characterizations from extreme surfaces, with the support of our GDMS instrumentations and more than 25 years of expertise in the field.

The complementarity of GDOES also relies on its capacity to analyze the evolution of H, C, N & O contents from the surface to the bulk of the material, qualitatively, and quantitatively in a second step.

ISO Standard Testing

GDOES is already known and used for standard testing, for example ISO 14707:2021 and ISO 16963:2017. The latter's field of application is presented below to show the wide range of contents that can be studied with GDOES.

- Zn from 0,01 wt% to 100 wt%;
- Al from 0,01 wt% to 100 wt%;
- Ni from 0,01 wt% to 20 wt%;
- Fe from 0,01 wt% to 20 wt%;
- Si from 0,01 wt% to 10 wt%;
- Pb from 0,005 wt% to 2 wt%;
- Sb from 0,005 wt% to 2 wt%.




Figure 3: Example of GDOES qualitative analysis to characterize the surface coating on a Steel sample

Qualitative and Quantitative Results

We have the capabilities to provide you with the most accurate analytical solution, while adopting a project-based approach with a first qualitative GDOES analysis, followed if necessary, by a quantitative strategy of external calibration with the appropriate Reference Materials.

Elemental Analysis by Glow Discharge Mass Spectrometry (GDMS)

Best Instruments for Best Quality

To stay at the highest technological level and to provide you with the best services, we seek to constantly develop our instrumentation and analytical services. Consequently, we now have eight instruments for GDMS services. These instruments allow us to accept a high number of samples, and to perform analyses quickly.




Figure 4: GDMS instruments

Different geometries of GDMS instruments are available within our French entity. This allows our specialists to continuously build expertise and know-how about their enhanced capabilities.




Figure 5: Simplified scheme of the GDMS principle

High Resolution GDMS

To keep it simple, High Resolution GDMS relies on the creation of a source of ions that are sequentially separated via Mass Spectrometry to accurately characterize sub-ppmwt contents in inorganic materials.

An Ultra-Pure Argon stream is introduced in the analytical cell where the sample has been first inserted. Voltage is applied between the anode and the sample in order to ignite the plasma. This Argon plasma bombards the surface of the sample. The sputtered species are then atomized, ionized, and extracted from the cell to fly through the High-Resolution (HR) analyzer system.

The HR analyzer is a double focused mass spectrometer (magnetic and electrostatic) which separates the ions depending on their mass on charge ratio m/z . The detection system allows the sequential determination of majors, minors, trace, and ultra-trace contents, down to the ppb level.

Sample Shape

HR-GDMS can be useful on all kinds of inorganic samples, conductive, semi-conductive and even non-conductive.

Solid conductive samples: the most straightforward way to analyze these samples is to introduce them directly in the instrument.

- As is, if the analytical surface is smooth and flat, and measures at least 20 x 20mm.
- Under pin form, if the sample can be cut at the following size: 20 x 2 x 2 mm. The analysis requires a preliminary mechanical preparation and a chemical etching to remove every potential contamination brought by the mechanical preparation.

Powders and non-conductive samples: The analysis of powders, or samples with a low conductivity, needs to be performed with the use of a high purity conductive binder (usually 7N grade Tantalum or Indium). This binder holds the sample mechanically and is also used as a secondary electrode to allow the ignition of the plasma.

Element	Concentration (ppmwt)	Element	Concentration (ppmwt)
Li	0.001	Ag	0.005
Be	0.001	Cd	0.01
B	0.005	In	0.005
F	0.05	Sn	0.005
Na	0.005	Sb	0.005
Mg	0.005	Te	0.05
Al	0.005	I	0.05
Si	0.001	Cs	0.005
P	0.005	Ba	0.005
S	0.05	La	0.005
Cl	0.01	Ce	0.005
K	0.01	Pr	0.005
Ca	0.01	Nd	0.005
Sc	0.001	Sm	0.005
Ti	0.005	Eu	0.005
V	0.001	Gd	0.005
Cr	0.005	Tb	0.005
Mn	0.005	Dy	0.005
Fe	0.005	Ho	0.005
Co	0.005	Er	0.005
Ni	0.005	Tm	0.005
Cu	Matrix	Yb	0.005
Zn	0.01	Lu	0.005
Ga	0.05	Hf	0.005
Ge	0.05	Ta	1
As	0.005	W	0.005
Se	0.01	Re	0.005
Br	0.005	Os	0.005
Rb	0.005	Ir	0.005
Sr	0.005	Pt	0.005
Y	0.005	Au	0.01
Zr	0.005	Hg	0.01
Nb	0.005	Tl	0.001
Mo	0.005	Pb	0.005
Ru	0.005	Bi	0.005
Rh	10	Th	0.0001
Pd	0.005	U	0.0001

Table 1: Typical reporting limits for a Copper sample analyzed by GDMS



ISO 9001 and 17025 accredited

Copyright © 2023 Eurofins Scientific | Rev: 1.03.2023 M-053921 EAG.COM

eurofins

EAG
Laboratories

Applications, Accuracy, and Uncertainty

HR-GDMS leads to the characterization of the elemental trace concentrations of the whole periodic table (Hydrogen and noble gases taken apart). High Resolution ensures the discrimination of most of the existing interferences, giving the ability to accurately measure percent to ultra-trace (ppbw) contents with a relative instrumental precision of about 20%.

COFRAC Accreditation - ISO 17025 Norm Certification for HR-GDMS Analyses

Our quality system and our continuous improvement policy has led our laboratory to be accredited by the French Accreditation Committee. Our HR-GDMS service has been accredited for most of the materials usually analyzed, for example steels, Nickel alloys, Titanium alloys, Aluminum alloys or Copper samples.

The accreditation scope is available on the COFRAC website (www.cofrac.fr) with the accreditation number 1-1993.




Figure 6: Simplified scheme of GDMS measurements

Interstitial Gas Analysis (IGA)

We are internally equipped with five IGA instruments, covering specific needs in the determination of interstitial gases H, N/O and C/ S, and completing the range of the results provided by GDMS and ICP-OES/MS. These instruments allow us to provide you with a complete evaluation of the whole periodic table, from Hydrogen to Uranium.




Figure 7: Simplified scheme of IGA measurements

Analysis of O, N & H Contents by Inert Gas Fusion

Inert gas fusion (using Argon or Helium) of a 1-gram sample is performed in a graphite crucible, deposited in an oven. The concentrations of the H₂, N₂ and O₂ gases that are released during the fusion are then determined using either a Thermal Conductivity Detector for H₂ and N₂, or an Infrared detector for O₂. At least three replicates per sample are analyzed, leading to a more representative evaluation of the whole sample.




Figure 8: LECO ONH836 Analyzer

Analysis of C and S Contents by Combustion

This time, combustion using pure Oxygen of a 1-gram sample is performed in a ceramic-made crucible. Concentrations of the gases released in the oven, CO₂ and SO₂, are then measured using an Infrared detector. At least three replicates per sample are also performed in order to have a more representative idea of the homogeneity of the sample



ISO 9001 and 17025 accredited

Copyright © 2023 Eurofins Scientific | Rev: 1.03.2023 M-053921 EAG.COM



Figure 9: LECO CS844 Analyzer

Applications and Limits of Detection (LOD)

All inorganic materials are theoretically suitable for IGA measurements. We can then consider the analysis of metallic matrices (Si, Cu, Ni, Co), alloys (Nickel superalloys), solid solutions (steels, shape memory alloys NiTi), refractory materials (SiC, Ta, W), oxides (Al₂O₃), nitrides (Si₃N₄), and many more.

Element	Ultimate Limit of Detection (ppmwt)
H	0.1
C	1
N	1
O	1
S	1

Table 2: Ultimate limits of detection obtained on 7N purity copper sample

ISO 17025 Accreditation for H, N, O, C and S Analyses

Our quality system and our continuous improvement policy has led our laboratory to be accredited by the French Accreditation Committee. Our IGA service received this accreditation for most commonly analyzed materials: copper, steels, nickel alloys, titanium alloys, and aluminum alloys. The accreditation scope is available on the COFRAC website (www.cofrac.fr) with the accreditation number 1-1993.

A High Standard Quality of Analysis

Our analytical methods are developed internally, while keeping American Society for Testing and Materials (ASTM) standards as models and examples.

- ASTM E-1447: H in Ti and Ti alloys
- ASTM E-1019: C, S, N & O in steels, Fe, Ni, Co and similar alloys

A daily quality control is performed with high metrological quality Reference Materials (CRM, primary standards) chosen from recognized providers like EURONORM-MRC, JSS, or BAM.

Element	Concentrations in the scope (ppmwt)
C	100 to 4000
N	30 to 300
O	10 to 150
S	5 to 250
H	0.5 to 10

Table 3: Elemental concentrations covered by our ISO 17025 accreditation scope



ISO 9001 and 17025 accredited

Copyright © 2023 Eurofins Scientific | Rev: 1.03.2023 M-053921 EAG.COM



EAG
Laboratories

Thermogravimetric and Differential Thermal Analyses (TGA/DTA)

Thermogravimetric Analysis, also known as TGA, monitors the evolution of a sample's mass as a function of temperature and time, under controlled atmosphere. TGA is very useful when precise information is needed concerning the determination of volatile elements or molecules, thermal stability, degradation characteristics and kinetics, or even behavior during sintering.




Figure 10: Simplified principle of a TGA/DTA measurement

For Differential Thermal Analysis, or DTA, the goal is to monitor the difference of temperature between the sample of interest and an inert reference, difference resulting from internal thermal events within the material. These internal events can be the result of a change of the chemical composition, or the crystalline structure of the studied material.

A new instrument has been acquired, Setaram LABSYS evo. It adds the possibility to monitor phenomenon from Room Temperature to 1600°C on organic, metallic or even ceramic materials. Atmosphere can be designed to fit your specific needs.




Figure 11: SETARAM LABSYS evo

All solid samples can be analyzed, without any preliminary preparation of the sample, including powders, chunks, flakes, and similar. The analysis can either be qualitative or quantitative.

We are ready to study your demand and adapt our existing procedures to fit your particular needs, to provide you with answers to your specific problems. This can include:

- Composition study, deformation
- Melting points, Lost on Ignition (LOI)
- Study of chemical reactions (hydrogenation; reaction kinetics; pyrolysis; sintering; ...)
- Monitoring of degradation, study of thermal stability
- Comparison between different batches, failure analysis, and so much more




Figure 12: Results obtained by TGA/DTA



ISO 9001 and 17025 accredited

Copyright © 2023 Eurofins Scientific | Rev: 1.03.2023 M-053921 EAG.COM

Analysis of Major Elements by Inductively Coupled Plasma Optical Emission Spectroscopy – ICP-OES

GDMS and IGA measurements are completed in our facility by Inductively Coupled Plasma Optical Emission Spectroscopy (ICP-OES), in order to accurately measure the major concentrations in the material of interest. ICPOES is a multi-element technique, relying on wet chemistry, as the samples need to be introduced in the plasma in the liquid form.




Figure 13: Simplified scheme of ICPOES measurements

We recently acquired a new instrument, Analytik Jena PQ9000, to perform these analyses.



Figure 14: Analytik Jena PQ9000 instrument

It is mandatory to first mechanically prepare the sample. Our laboratory is equipped with instruments and tools designed to mechanically prepare all types of samples. It is then possible to transform solid samples into small chunks, or sample special areas of interest. The samples are then digested using specific acid mixtures, that are internally developed as a function of the material needed to be analyzed.

Once again, specific equipment, such as hot plates, microwaves or even high-pressure bombs in an oven, is used to digest the most difficult matrices.

Instrument Check and External Calibration Curves

Analysis of the samples is performed after a daily instrumental check, against a list of criteria, as for example signal intensity, stability and oxides level.

The instrument is then calibrated to measure its response as a function of the analyte's concentration. Elemental standard solutions, traceable to NIST, are used to gravimetrically prepare calibration curves, allowing an increased accuracy.




Figure 15: Microwave-assisted digestion system (left) and instrumental calibration curves (right)

Applications and Strengths

Our ICPOES instrument is dedicated to the analysis of major elements in organic or inorganic materials. Our laboratory is able to deal with metallic alloys, superalloys, oxide powders, and many more. Our experienced team can adapt a solution to your particular sample, if your sample request more than a straightforward digestion or analysis.

ICPOES enables the determination of concentrations associated to instrumental uncertainties of about 5% for concentrations above 1 wt%, and about 10% for concentrations below 1 wt%.

ISO 9001 and 17025 accredited



Copyright © 2023 Eurofins Scientific | Rev: 1.03.2023 M-053921 EAG.COM

 eurofins

EAG
Laboratories

Our internal procedures include the sampling of 3 different replicates (generally between 0.1 and 0.5g) from the initial sample. The average of the 3 measurements is reported, with the standard deviation of these replicates. This information can be used to understand the homogeneity of the sample if the standard deviation is low, or inhomogeneity if the standard deviation is higher.

Element	Concentration (in wt%)	Specifications (in wt%)
Si	0.381+/-0.006	< 0.75
Cr	19.50+/-0.08	18-20
Mn	1.76+/-0.01	2 max
Fe	Balance	Balance
Ni	9.71+/-0.03	8-12

Figure 17: Example of a steel sample: ICPOES results and specifications

The example presented in Figure 17 shows that the ICPOES measurements are fulfilling client's expectations, even for the tightest specifications. A closer look on the standard deviations show a low dispersion of the results coming from the 3 replicates, traducing the homogeneity of the analyzed sample.

Contact Us

For more information concerning the techniques presented in this document, or the complementary techniques provided by our Eindhoven, Netherlands location or any of our US sites, please feel free to contact us, using the following phone number and email address. We'll be pleased to find a dedicated solution for you.

+33 (0)5 61 73 15 29

info.fr@eag.com

Accréditation n°1-1993
Portée d'accréditation disponible sur www.cofrac.fr



ISO 9001 and 17025 accredited

Copyright © 2023 Eurofins Scientific | Rev: 1.03.2023 M-053921 EAG.COM