

IMPROVED PROFILE ACCURACY USING **POINT-BY-POINT CORRECTED-SIMS**. A TECHNOLOGY EXCLUSIVE TO EAG.



EAG
Laboratories



What is PCOR-SIMSSM?

PCOR-SIMSSM is a technique that can measure layer thickness, composition, and doping profiles more accurately than a regular SIMS analysis. It uses the simple idea that composition in a multi-layer system can be determined by the SIMS measurement itself. If composition can be determined then the correct sputter rates and sensitivity factors for impurities can be applied at every point.

Where did PCOR-SIMSSM Originate From?

It was commonly assumed that Secondary Ion Mass Spectrometry (SIMS) could not quantify matrix-level concentrations with sufficient accuracy, and calibrating dopant intensities continuously, with respect to varying matrix compositions, as it would be prohibitively time consuming. Scientists at Eurofins EAG Laboratories recognized these challenges, and the need for a deeper level of accuracy, so they perfected the methodology and technique to birth PCOR-SIMSSM, or Point-by-Point CORrected-SIMS.

How is Composition Measured?

Standards need to be prepared covering the range of expected composition for a given material system. From profiles on these standards, calibration functions for composition and sputter rate can be determined.

Calibration Functions for Impurities

Known amounts of impurities are introduced into a set of composition standards. Sensitivity factors for impurities can then be determined as a function of composition.



What is SIMS?


SIMS is a powerful depth profiling tool with unmatched sensitivity to all elements. However if the composition of the system changes with depth, then concentration profiles become less accurate due to 'matrix effects'!

What are Matrix Effects?


In SIMS, a primary ion beam sputter etches a surface. Only a fraction of the sputtered material is ionized during sputtering. The ion fraction depends on the element and also on the surface chemistry - the matrix material. Change the matrix material and the ion yields change.

Why are Matrix Effects Important?

If composition changes with depth, then dopant and contaminant ion yield can change. Concentration profiles may not be accurate. Sputter rate can also change resulting in a less accurate depth scale.



Sputter Rates in Silicon Based Materials
Sputter rates vary by a factor of 1.7. Data are normalized to 1 for silicon.




SIMS Matrix Effect

Ion yield for silicon is lower in Si than in SiO₂. The profile shows Si ion counts decrease from SiO₂ to Si even though concentration is actually 3x higher!

The boron profile shows a drop-out in the Si layer but the actual boron implant is continuous as indicated by the dotted line.




Applications of PCOR-SIMSSM for Multi-Layer Systems




VCSEL n-DBR Layers

PCOR-SIMSSM can determine the layer composition at every data point. In this case we show a GaAs/AlGaAs multi-layer structure. We can then plot a silicon dopant profile using the correct silicon sensitivity factor for the composition at every data point (shown in green). Silicon plotted using a single GaAs sensitivity factor is shown in red.



pHEMT


Many pHEMT structures benefit from PCOR-SIMSSM depth profiling. Composition for the GaAs, AlGaAs, and InGaAs layer structure is determined at every data point. Accurate layer thickness is shown at every data point since sputter rate variations for each composition is taken into account. Silicon spike concentration and spacer thickness can be accurately measured. Comparison with TEM shows the alignment of the dopant spikes with the sample structure.



Ultra-Shallow Implants

Very shallow-high dose implantation results in a significant portion of the dopant residing with a surface modified layer. This dopant is not available for activation.

PCOR-SIMSSM shows the thickness of the surface modified layer, measures oxygen content, and determines the correct sensitivity factor at each data point for the dopant. The result is an accurate dopant profile shape and dose, and a more realistic measure of junction depth, defined by the distance from the oxide/Si interface to the point where the boron profile crosses 5E18 at/cc.




SiGe

PCOR-SIMSSM plotted in green shows the correct SiGe composition profile, the correct carbon dopant and oxygen contaminant profiles, and the correct thickness at every data point. The red plot shows uncorrected composition, carbon and layer thickness where carbon and oxygen are 2x too high and layers are 20% too thin.

How Do We Know PCOR-SIMSSM is Correct?

At EAG we continuously validate PCOR-SIMSSM against reference materials calibrated using other techniques.





GaAs/AlGaAs Layers

Measurement of composition is the first step in PCOR-SIMSSM. Here we compared PCOR-SIMSSM to RBS (Rutherford Backscattering Spectrometry), a technique that is not matrix dependent and does not need standards. The results agree to within 5%.

B High Dose Implant Through Surface Oxide

PCOR-SIMSSM is needed to extend accurate dopant quantification into percent level. Here we compare to ERDA (elastic recoil detection analysis) with PCOR-SIMSSM. ERDA provides a quantitative depth profile in the high concentration portion of the profile that is not matrix dependent and does not need standards. ERDA confirms the PCOR-SIMSSM profile is correct.

ERDA data courtesy of Dr. Wilfried Vandervorst.



B Plasma Implant into Polysilicon

PCOR-SIMSSM is needed to accurately measure dopants in multiple layers and blended composition layers. Here we compare PCOR-SIMSSM with TEM and NRA (nuclear reaction analysis). TEM confirms the location and thickness of layers and NRA confirms that total boron is correct.



About Eurofins EAG Laboratories

When it comes to understanding the physical structure, performance, chemical properties and composition of materials, no other scientific services company offers the breadth of experience, diversity of analytical techniques or technical ingenuity of Eurofins EAG Laboratories. We don't just perform testing, we drive commercial success—through thoughtfully designed investigations, technically superior analyses, and expert interpretation of data.

We deliver multi-disciplinary, problem-solving expertise to help our customers accelerate innovation, ensure quality and safety, and protect intellectual property. Whether you are seeking to reduce time-to-market, solve manufacturing problems or ensure regulatory compliance, turn to Eurofins EAG. We know how to bring the power of science to every phase of your product lifecycle.

- 20+ facilities located in the US, Europe, and Asia
- 2,500+ instruments
- 1,000+ highly-educated employees
- Serving more than 5,000 clients worldwide
- Revenue sourced from more than 50 countries

+1 800-366-3867

WWW.EAG.COM

Complete Lifecycle of Services

 PRODUCT INNOVATION & IMPROVEMENT

 INVESTIGATION & TROUBLESHOOTING

 QUALITY ASSURANCE

 MANUFACTURING SUPPORT

 FAILURE ANALYSIS

 REGULATORY COMPLIANCE

 MANUFACTURING & SUPPLY CHAIN SUPPORT

 CONSULTING & LITIGATION

 eurofins

EAG
Laboratories